

Application of a real-time polymerase chain reaction assay for the myxozoan parasite *Henneguya ictaluri* and its year-round prevalence in catfish aquaculture

\*Matt J. Griffin, David J. Wise, Linda M. Pote, Alvin C. Camus, Michael J. Mauel, Terrence E. Greenway

Thad Cochran National Warmwater Aquaculture Center  
Mississippi State University  
Stoneville, MS 38776



Author/ISAAH

*Henneguya ictaluri*


- Myxozoan parasite
- Causative agent of proliferative gill disease (PGD) in channel catfish
  - Hamburger gill
- 3<sup>rd</sup> most commonly diagnosed disease at the Aquatic Diagnostic Lab, Stoneville, MS
- Mortalities can approach 100% in severe outbreaks



Author/ISAAH

© Author/ISA AH  
*Henneguya ictaluri*

- Life cycle involves ubiquitous benthic oligochaete *Dero digitata*
- Exposure to the actinospore stage results in a severe inflammatory response in channel catfish gills causing respiratory distress
- Currently there are no treatments
- Primarily a problem in the spring, with lesser incidence in the fall, when water temps are between 16-25 °C



Bar = 20 µm

3

© Author/ISA AH  
*Henneguya ictaluri*

- Mostly an issue for fingerlings stocked newly into ponds for grow-out
- Can cause significant losses larger fish, but most severe outbreaks occur in smaller, younger fish
- Different populations of fish in the pond can be affected differently
- Effects age, size and previous exposure to the organism have on outbreak severity is poorly understood
  - Fish appear to be able to acquire resistance
- Once fish are removed from the source of infection they quickly recover

4

## © Author/ISAAH Problem

- Following an outbreak which results in significant losses, when is it safe to re-stock fish into the pond?
  - Resident fish often show no signs of disease yet pond actinospore concentrations are lethal to newly introduced fish
- During an outbreak, if moving fish is a feasible option how do we identify parasite free environments for which to relocate fish?
- Called for a method to identify approximate actinospore concentrations in a pond prior to stocking

## © Author/ISAAH Current practices

- Fish Health Program at NWAC utilizes sentinel fish (Wise et al. 2004; 2008)
  - Two separate week long exposures
  - Quantify damage by calculating the percent filaments exhibiting at least one chondrocytic lesion or “break”
  - Does it get better or worse between exposures?
  - Helps determine dynamics of a given outbreak.
  - Outbreaks do not tend to relapse
    - Once an outbreak is done, its done for that season

6

## © Author/ISAAH Disadvantages . . .

- Takes a minimum of one week to determine status of pond
- Takes a minimum of two weeks to determine the dynamics of the outbreak
- Requires equipment, source of PGD free fish, and is fairly labor intensive
- Is inconclusive if cage system fails
  - Lose cage/fish
  - All fish die before sampling

7

## © Author/ISAAH Alternatives

- Using QPCR can we measure the actinospores directly from pond water samples?
  - Eliminates need for sentinel fish
  - Significantly less labor intensive than sentinel fish exposures
  - Much more rapid than sentinal fish exposures

8

How do you get from this . . .



to this?



9



10

## Water sampling, DNA isolation and PCR


- Sorvall RC-6 Ultracentrifuge
  - A series of centrifugation steps concentrates >500 mls of pond water and sediment into 1.5 ml tubes
- Mo Bio Soil DNA Isolation Kits
  - Powersoil
  - Ultraclean
- Bovine Serum Albumin added to PCR reaction
  - 1 µg/µl final reaction concentration
- ABI's Environmental Taqman Mastermix 2.0

11

## Pond trials

- 20 replicate pond trials
- Sentinel fish were held in pond for seven days
- Water was sampled from two different geographic locations within the pond in coordination with the placement and collection of sentinel fish
- % mortality was correlated to actinospore equivalents per liter of pond water

12



## Conclusion . . .

The *H. ictaluri* specific QPCR assay can detect *H. ictaluri* in environmental samples and has tremendous potential in both research and production settings.


The ability to quantify *H. ictaluri* DNA from fish tissue and environmental samples provides an invaluable tool in diagnostics, epidemiological studies, treatment evaluations and management of this disease.

## Year-round prevalence of *H. ictaluri*

- In the spring of 2009 a producer was experiencing significant PGD related mortalities on his operation
- Once monthly we sampled pond water from 24 ponds on his operation
  - (May 2009-May 2010)

15


## Percent of PGD cases submitted by month (2007-2009)




16



## % Positive Ponds (by month)



## Mean actinospore equivalents from positive ponds



## © Author/ISAAH Conclusions

- *H. ictaluri* is present in catfish ponds year round, but does not always cause disease
- *H. ictaluri* is present in numbers that would result in significant mortality throughout the year and not only in the spring
- Previous exposure to the organism may provide protection against re-exposure

19

## © Author/ISAAH New questions . . . .

- What effect does age, size, and or previous exposure have on PGD severity
  - Takes a larger dose of parasites to affect larger fish, which are normally older, and in the pond, more likely to have seen *H. ictaluri* before
- What effect does temperature have on PGD?
  - If PGD is actually an immune response to *H. ictaluri*, in an immune suppressed state (cold temperatures) do fish get PGD?
- If previous exposure during an immune suppressed state provides protection, would fall stocking reduce incidence of PGD in the spring?

20